

STUDY REPORT

MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS AND OTHER NON-POROUS SURFACES

ISO 22196:2011

Page 1 of 6

Job reference — J002048-2

Lab Ref/Report No.	J002048
Testing Laboratory Site	Microbiological Solutions Limited Gollinrod, Bury, BL9 5NB
Company owner	Angela Davies, Managing director
Report Date	04/08/2020
Period of Analysis	23/07/2020-25/07/2020

Customer	Hale Environments Ltd
Contact Name	Giles Tickle
Address	James House, Mere Park, Dedmere Road, Marlow, SL7 1FJ
Email	info@hale-environments.co.uk
PO Number	Q003134

Name of product	Hale Hygiene +
Batch number	N/S
Storage Conditions	Ambient
Appearance of the Product	White paint
Method	ISO 22196:2011
Neutraliser	N1
Volume of test inoculum used	0.4ml
Dimensions of test specimens	76mmx26mm
Dimensions of cover film	66mmx16mm square thickness 0.065mm
Cover film material	Polythene
Test Temperature	20°C
Temperature of Incubation	Bacteria - 37°C ±1°C for 24hr to 48hrs
Identification of the reference strains	<i>Escherichia coli</i> NCTC 10418 (ATCC 10536) <i>MRSA</i> NCTC 12493

STUDY REPORT

MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS AND OTHER NON-POROUS SURFACES

ISO 22196:2011

Page 2 of 6

Job reference — J002048-2

Introduction

ISO 22196:2011 specifies a method of evaluating the antibacterial activity of antibacterial-treated plastics, and other non-porous, surfaces of products (including intermediate products).

It is not intended to be used to evaluate the effects and propagation of bacteria on non-porous surfaces without antibacterial treatments. ISO 846 describes tests to evaluate the effects and propagation of bacteria on non-porous surfaces, which are different from those covered by ISO 22196.

Secondary effects of antibacterial treatments, such as the prevention of bio deterioration and odour, are not covered by the standard, which is not intended to be used or referenced as a method to document or claim biodegradability of, for instance, plastics materials.

Building materials are excluded, except where they are used in the same manner as treated articles.

Antibacterial-treated textile products are excluded, even if the surfaces are coated or laminated (such products are covered by ISO 20743).

Photocatalytic materials and products are excluded (such materials and products are covered by ISO 27447).

Test Method – Standard conditions

The test specimen is placed in an empty agar dish and 0.4ml of the test inoculum pipetted onto the test surface. The inoculum is then covered with a 40mm X 40mm square piece of film and gently pressed down ensuring the liquid does not leak beyond the edge of the film. The surface is maintained at $35^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 24 hours ± 1 hr. Testing is performed in triplicate with organisms recovered from control samples at 0 hours and 24 hours and from test samples at 24 hours only. Organisms are recovered by washing test specimens in 10ml of validated neutralizer and enumerating the wash liquid by pour plate.

Deviations from standard method

The product was tested against E.coli and non-standard organism MRSA.

The product was applied to glass slides prior to testing. The size of the glass slides was factored into final calculations.

Plain glass slides were used as controls.

STUDY REPORT

MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS AND OTHER NON-POROUS SURFACES

ISO 22196:2011

Page 3 of 6

Job reference — J002048-2

Acceptance Criteria

The value of the antibacterial activity can be used to characterize the effectiveness of an antibacterial agent. The antibacterial-activity values used to define the effectiveness shall be agreed upon by all interested parties.

Conclusion

The test product has shown a 3.14 (99.928%) log reduction against E.coli and a 3.49 (99.968%) log reduction against MRSA, when tested according to the conditions stipulated in this report.

Test results were calculated against a plain glass control material without paint.

See raw data tables below for test results.

The sample will be retained for 1 month unless otherwise requested.



Laboratory Manager
Megan Barrett



Technical Project Manager
Peter Thistlethwaite

The test results on this report refer only to the items tested. This report shall not be reproduced except in full and with written approval of Microbiological Solutions Ltd

STUDY REPORT
MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS
AND OTHER NON-POROUS SURFACES
ISO 22196:2011

Job reference – J002048-2

Test results – MRSA

24 hours					
MRSA	Test	Dilution	cfu/ml	Average cfu/ml	At (Log10 bacteria/cm2)
		0	8	10	0.98
		0	19		
		0	3		

0 hour					
MRSA	Control	Dilution	cfu/ml	Average cfu/ml	U0 (Log10 bacteria/cm2)
		3	32	48	4.65
		3	21		
		2	90		

24 hour						
MRSA	Control	Dilution	cfu/ml	Average cfu/ml	Ut (Log10 bacteria/cm2)	Antibacterial activity (R)
		3	33	31	4.46	3.49
		3	42			
		3	17			

STUDY REPORT
MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS
AND OTHER NON-POROUS SURFACES
ISO 22196:2011

Job reference – J002048-2

Test results – E.coli

24 hour					
Escherichia coli	Test	Dilution	cfu/ml	Average cfu/ml	At (Log10 bacteria/cm2)
		0	38	23	1.33
		0	12		
		0	18		

0 hour					
Escherichia coli	Control	Dilution	cfu/ml	Average cfu/ml	U0 (Log10 bacteria/cm2)
		2	92	81	3.89
		2	79		
		2	73		

24 hour					
Escherichia coli	Control	Dilution	cfu/ml	Average cfu/ml	Ut (Log10 bacteria/cm2)
		3	32	31	4.47
		3	42		
		3	19		

Antibacterial activity (R)
3.14

STUDY REPORT
MEASUREMENT OF ANTIBACTERIAL ACTIVITY ON PLASTICS
AND OTHER NON-POROUS SURFACES
ISO 22196:2011



Page 6 of 6

Job reference – J002048-2

Key

I: - Test inoculum concentration cfu/ml

U₀ – Log10 cfu/cm² recovered from control at time point 0h

A_t – Log10 cfu/cm² recovered from test sample at time point 24h

U_t – Log10 cfu/cm² recovered from control at time point 24h

R = (U_t – U₀) – (A_t – U₀) = U_t – A_t